



Glycohemoglobin Analyzer HbA1c BioHermes A1C EZ 2.0

Easy and accurate portable testing of long-term blood glucose level HbA1c, for better diabetes monitoring in hospitals and at home.

Glycohemoglobin Analyzer HbA1c

Convenient

- Only one drop (~ 3 microliter) of fingertip blood needed
- Easy to operate in three steps
- Voice prompt during the measuring process
- Handheld, battery-driven, portable
- Automatic storage of 1,000 test results with corresponding test time and date
- Mini USB interface for easy data transfer

Fast

- No calibration required
- Test results within 5 minutes
- No preheating needed before testing

Precise

- Accurate results with CV <3%
- NGSP and IFCC double certificates, CE certified

Results displayed as:

- % HbA1c (DCCT - preset) or mmol/mol (IFCC)
- Additional display of eAG in mmol/L

Accessories:

- 4 AAA Batteries
- User manual, operation guide, cleaning and maintenance guide

Kit de test:

- 25 Test strips
- Buffer A | Buffer B
- 25 Blood sampler
- Code chip (for matching the analyzer with the test strips)

Test Kit



Immediate feedback
in only 5 minutes



For more information, please contact us. We remain at your service!

action medeor e.V.
German Medical Aid Organization

St. Töniser Str. 21
47918 Tönisvorst | Germany
Phone: +49 2156 9788 100

info@medeor.de | www.medeor.de/en

* Please note that shipping including batteries is not allowed per postal services.